

FDM ATR Retail Adhesives & Sealants

Epoxy Resin



Epoxy Hardener



Epoxy Cured



Chemical Class	Applications:
Acrylates	Arts & Crafts
Butyl Rubbers	Asphalt
Cyanoacrylates	Automotive
Elastomers	Construction
Expoxies	Contact Cement
Intumescent	Gasket Sealants
Polyethers	Kitchen & Bath
Silicate Cements	Marine
Siliconess	Masonry
Siliconized Caulks	Rubber Cement
Urethanes	Tacky Glue
Urethanized Acrylics	Threadlockers
Vinyls	Wood Glue
...and more	...and more

Features:

- 458 spectra
- Monolithic diamond ATR crystal
- 4000 to 400 cm⁻¹
- Complements the FDM ATR Polymers and FDM Raman Retail Adhesives and Sealants

Samples:

3M Auto Plastic Emblem and Trim Adhesive
 3M Carpenters Wood Glue, Indoor
 3M Ceramic Adhesive, Outdoor, cured...

ELMERS® China plus Glass Cement
 ELMERS® Contact Cement
 ELMERS® Glue-All® Max...

Gold Eagle® AlumAseal® Radiator Stop Leak
 Goo Gone® Pro Power™, PNGG12DO
 Gorilla Glue, Precision Glue Pen



Fiveash Data Management, Inc.

608-236-9145

www.fdmspectra.com

Fax: 608-236-9170