



FDM Raman Polymers

- Polymers
- Plasticizers
- Monomers
- Trade name polymer additives
- Common and engineering plastics
- Chemical and industrial names
- Many CAS numbers
- Many Mw, Mn and Tg values
- Chemical names with CAS numbers
- 780 nm laser, 3250 cm⁻¹ to 200 cm⁻¹
- White Light Intensity Correction
- Complements the FDM ATR Polymers and FDM Raman Adhesives & Sealants.


Tripropylene glycol diacrylate,
42978-66-5




1,2-Polybutadiene,
Mw=100000, 31567-90-5



Dicyclohexyl phthalate,
84-61-7



Acrylonitrile butadiene copolymer,
butadiene 45 percent, 9003-56-9



Fiveash Data Management, Inc.

608-236-9145

www.fdmspectra.com

Fax: 608-236-9170